MONTH &	YEAR:	

San Diego Miramar College Automated External Defibrillator AED CHECKLIST

Serial Number:	
Building:	

DAILY CHECK	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
Unit is present with ancillary pouch & spare battery																															
Check status indicator on AED handle for "OK"	1																														
Verify daily check with INITIALS																															

MONTHLY CHECK

Examine AED case, connector and battery																
Examine battery pins for bending or discoloration																
Check expiration date on battery and all electrodes																
With battery installed: Turn on AED																\Box
Self Test Message																
Momentary illumination on all LED and all LCD segments																
Battery low or replace																
Battery self-test																
Service indicator-call service message	l I															
Check contents of ancillary pouch for razor, scissors, towel and pocket mask																
INITIAL Monthly-Check Complete																